

2

DTIC FILE COPY

MEMORANDUM REPORT BRL-MR-3866

BRL

AD-A228 213

IMPROVING SPECTRAL FITS OF ABSORPTION DATA
TAKEN WITH AN ARRAY DETECTOR:
WAVELENGTH "LINEARIZATION"

JOHN A. VANDERHOFF
ANTHONY J. KOTLAR

SEPTEMBER 1990

DTIC
ELECTE
OCT 29 1990
S E D

APPROVED FOR PUBLIC RELEASE; DISTRIBUTION UNLIMITED.

U.S. ARMY LABORATORY COMMAND

BALLISTIC RESEARCH LABORATORY
ABERDEEN PROVING GROUND, MARYLAND

90 10 24 032

NOTICES

Destroy this report when it is no longer needed. DO NOT return it to the originator.

Additional copies of this report may be obtained from the National Technical Information Service, U.S. Department of Commerce, 5285 Port Royal Road, Springfield, VA 22161.

The findings of this report are not to be construed as an official Department of the Army position, unless so designated by other authorized documents.

The use of trade names or manufacturers' names in this report does not constitute indorsement of any commercial product.

UNCLASSIFIED**REPORT DOCUMENTATION PAGE**Form Approved
OMB No. 0704-0188

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing the burden, to Washington Headquarters Service, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302, and to the Office of Management and Budget, Paperwork Reduction Project (0704-0188), Washington, DC 20503.

1. AGENCY USE ONLY (Leave blank)		2. REPORT DATE September 1990	3. REPORT TYPE AND DATES COVERED Final, Mar 90 - Aug 90	
4. TITLE AND SUBTITLE IMPROVING SPECTRAL FITS OF ABSORPTION DATA TAKEN WITH AN ARRAY DETECTOR: WAVELENGTH "LINEARIZATION"			5. FUNDING NUMBERS 1L161102AH43	
6. AUTHOR(S) John A. Vanderhoff Anthony J. Kotlar			8. PERFORMING ORGANIZATION REPORT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES)				
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(ES) Ballistic Research Laboratory ATTN: SLCBR-DD-T Aberdeen Proving Ground, MD 21005-5066			10. SPONSORING / MONITORING AGENCY REPORT NUMBER BRL-MR-3866	
11. SUPPLEMENTARY NOTES				
12a. DISTRIBUTION / AVAILABILITY STATEMENT Approved for public release; distribution unlimited			12b. DISTRIBUTION CODE	
13. ABSTRACT (Maximum 200 words) Wavelength "linearization" is incorporated into a least squares fitting routine tailored for optical absorption spectra taken with a photodiode array. In general, the wavelength linearization algorithm improves the quality of the fits; in some cases, the standard deviation of the fits is improved by more than a factor of two.				
14. SUBJECT TERMS Absorption, Solid Propellant, Combustion, Temperature, Flame, OH Concentration, Photodiode Array, Monochromator, Wavelength Linearization. J50			15. NUMBER OF PAGES 22	
			16. PRICE CODE	
17. SECURITY CLASSIFICATION OF REPORT Unclassified	18. SECURITY CLASSIFICATION OF THIS PAGE Unclassified	19. SECURITY CLASSIFICATION OF ABSTRACT Unclassified	20. LIMITATION OF ABSTRACT UL	

NSN 7540-01-280-5500

UNCLASSIFIEDStandard Form 298 (Rev 2-89)
Prescribed by ANSI Std Z39-18
298-102

INTENTIONALLY LEFT BLANK.

TABLE OF CONTENTS

PAGE

LIST OF FIGURES	v
I. INTRODUCTION	1
II. EXPERIMENTAL	1
III. DATA ANALYSIS	2
IV. RESULTS	4
V. CONCLUSION	9
REFERENCES	11
DISTRIBUTION LIST	13

Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input checked="" type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By _____	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A-1	

INTENTIONALLY LEFT BLANK.

LIST OF FIGURES

<u>FIGURE</u>	<u>PAGE</u>
1. The light path for optical absorption measurements of OH	1
2. Configurations of the monochromator focal plane	2
3. A rotationally resolved OH absorption spectrum with the monochromator operated in second order and no linearization	5
4. Same conditions as Fig. 3 except here the monochromator was operated in first order	5
5. Same conditions as Fig. 3 except the data were fitted using wavelength linearization	6
6. Plot of the absolute difference of the observed values and the fitted values versus wavelength for the conditions of Fig. 3	7
7. Plot of the absolute difference of the observed and the fitted values versus wavelength for the conditions of Fig. 5	7
8. A rotationally resolved OH spectrum taken from 0.4 to 0.9 mm above the surface of a nitramine propellant burning in 1.5 MPa nitrogen; the monochromator was operated second order and the data were fitted using wavelength linearization	8

INTENTIONALLY LEFT BLANK.

I. INTRODUCTION

We have been conducting optical diagnostic measurements on solid propellant combustion for the past few years. While analyzing some recent OH absorption spectra^{2,3} with a least squares fitting program, an inconsistency appeared, *i.e.*, the computer fits of OH absorption data taken under high spectral resolution resulted in a larger statistical uncertainty than the fits for spectral data taken with half the spectral resolution. The source of this inconsistency has been identified and corrected using a wavelength "linearization"¹ method for an array detector.

II. EXPERIMENTAL

Premixed $\text{CH}_4/\text{N}_2\text{O}$ laminar flames and nitramine propellant flames have been probed by detecting rotationally resolved absorptions of the OH molecule around 306.4 nm. Details of the experimental apparatus appear elsewhere^{2,5} and only the optical detection segment will be described here. The optical path for light produced by a 1 kw mercury-xenon arc lamp is shown on Fig. 1. This broadband light is focussed and




Figure 1 The light path for the optical absorption experiments. The glass filter (GF) was used to eliminate first order light when operating the spectrometer second order; the sample is either a premixed flame or propellant.

apertured with lenses and adjustable irises. A 150 micron pinhole provides the final spatial restriction before reaching the sample region. Light transiting the sample region is line focussed onto a monochromator with a cylindrical lens. This Czerny-Turner type monochromator (model HR-320 manufactured by J-Y Optical Systems) has a 0.32 m focal length and a 2400 grooves/mm holographic grating. The monochromator was operated in second order to adequately resolve the rotational absorption lines occurring in the $\text{A}^2\Sigma^+ - \text{X}^2\Pi_1$ (0,0) vibrational band system of OH. First order light is rejected by a glass filter

(Schott UG-11) placed in front of the monochromator. A 0.025 mm entrance slit gives a spectral resolution (FWHM) of 0.03 nm for second order operation and a 0.06 nm resolution for first order operation.




Figure 2 Configuration where monochromator focal plane is (a.) normal to the central wavelength of and (b.) inclined at an angle G with respect to the central wavelength of the photodiode array.

The light detector is a photodiode array with 1024 photodiodes of 0.025 mm width. About 700 of these photodiodes receive amplified light from a microchannel plate intensifier. The geometry of the detector with respect to the diffraction grating of the monochromator (see Fig. 2) determines the non-linearity produced in the conversion of photodiode position to wavelength. The tilt angle (G) determines the extent of this non-linearity. In the data analysis section which follows, an algorithm is given for linearizing the wavelength.

III. DATA ANALYSIS

When a monochromator is used with an array detector, the relationship between wavelength and pixel (photodiode) position may not be linear. The following abbreviated data analysis describes such a case for a Czerny-Turner type configuration and is taken from an optical spectroscopy tutorial by Lerner and Thevenson.⁹

The grating equation can be written as

$$\sin A + \sin B = 10^{-6} knw \quad (1)$$

and

$$D_v = B - A \quad (2)$$

All of the symbols are shown in Figs. 2a and 2b and described in Table 1. The angle of

Table 1. Symbols, definitions, units and values when appropriate.

<u>Symbol</u>	<u>Definition</u>	<u>Units</u>	<u>Values for HR-320</u>
A	Angle of incidence	degrees	
B	Angle of diffraction	degrees	
D_v	Deviation angle	degrees	24
G	Tilt angle	degrees	2.4
k	Diffraction order	integer	
n	Grating groove density	grooves/mm	2400
w	Wavelength	nanometers	
w_c	Wavelength at center of array	nanometers	
w_m	Wavelength at position m	nanometers	
Bw_m	Angle of diffraction at wavelength m	degrees	
L_c	Focal length of monochromator	mm	320
L_{H1}	Perpendicular distance from grating to the focal plane	mm	319.719
B_{H1}	Angle from L_{H1} to the normal to the grating	degrees	
HBw_m	Distance from the intercept of the normal to the focal plane to the wavelength w_m	mm	
HBw_c	Distance from the intercept of the normal to the focal plane to the wavelength w_c	mm	13.4
Pwd	Nominal width of the pixel	mm	0.0245
Pm	Pixel number at wavelength w_m	integer	
Pc	Pixel number at wavelength w_c	integer	

incidence can be expressed as

$$A = \sin^{-1} \{ 10^{-6} knw / 2 \cos(D_v / 2) \} - D_v / 2. \quad (3)$$

For a specific wavelength, w_m , Eqn. (1) can be rewritten as

$$w_m = [\sin A + \sin Bw_m] 10^{-6} / kn \quad (4)$$

where

$$Bw_m = B_H - \tan^{-1}(HBw_m/L_H) \quad (5)$$

and

$$HBw_m = HBw_c - Pwd(Pm - Pc). \quad (6)$$

Eqn. (4) is now used to compute the wavelength for the least squares fitted OH spectra presented in the results section; previous published results^{2,5} used a fixed wavelength increment per pixel. The pixel width has been allowed to vary in the fitting procedure; although the nominal width between pixels is given as 0.025 mm the best fits occurred for values of about 0.0245 mm.

The main purpose of this diagnostic technique is to extract temperature and concentration from absorption spectra by using a non-linear least squares fitting procedure with an equation of the form

$$I_{tr} = \int S(w, w_o) I(w) dw \quad (7)$$

where I_{tr} , the integrated light transmitted, is a function of w_o , the index frequency of the diode; $S(w, w_o)$ is an instrument function appropriate for the spectral resolution of the detection system; and $I(w)$ is the transmitted intensity of a group of lines for a given species. A much more complete explanation of the absorption technique and the fitting procedure is described elsewhere.^{2,3,7,8}

IV. RESULTS

Figures 3 and 4 show spectral absorption data taken in an atmospheric pressure CH_4/N_2O laminar flame. These absorptions are from rotationally resolved transitions in the A-X (0,0) vibrational band system of OH. The data are represented by dots and the non-linear least squares fit to the data is a solid line. Conditions under which these absorption spectra were taken plus the fitted values obtained for the temperature and concentration are given in Table 2. The main difference† between the spectral data of Fig. 3 and Fig. 4 is that the monochromator was operating in second order for Fig. 3, therefore, these data are much better resolved. Contrary to intuition, however, the standard deviation of the overall fit and of both the temperature and concentration is higher for the better resolved

†It is also seen that the values of the temperature and concentration are different. This result is not surprising since there was no quantitative control on flame conditions and the data for Figs. 3 and 4 were obtained on different days.




Figure 3 Rotationally resolved, second order (uncorrected) OH absorption spectrum taken 1 mm above the burner surface of the premixed flame; solid line, least squares fit for values in Table 2.




Figure 4 Same conditions as Fig. 3 except here the monochromator was operated in first order.

data (see graph captions). This observation prompted our use of the wavelength correction algorithm discussed in the data analysis section. The data of Fig. 3 were subsequently refitted incorporating this algorithm and are plotted along with the wavelength corrected fit on Fig. 5. Although the standard deviation of the fit has dropped by more than a factor of two, it is difficult to see the differences in the quality of the fits when comparing Figs. 3




Figure 5 Same conditions as Fig. 3 except the data were fitted using wavelength linearization.

and 5. However, plots of the absolute difference (data point value - calculated value) for both Figs. 3 and 5 clearly show the improvement attained when incorporating this wavelength linearization algorithm, and these are given on Figs. 6 and 7, respectively. The absolute value differences shown in Fig. 6 (uncorrected) are about twice as large as in Fig. 7 (corrected). Moreover, the data of Fig. 7 appear random whereas the data of Fig. 6 are structured indicating the presence of a systematic error.

The data of Fig. 4 were obtained in first order for the same wavelength range covered on Fig. 3. This means that for the data of Fig. 4 the spatial extent of the dispersed light is about half that for the conditions of Fig. 3, hence, the nonlinearity should be less influential. This is confirmed by comparing the standard deviation of the least squares fit of the low resolution spectrum with and without the wavelength linearization. For spectra obtained with the grating operating in first order, no appreciable improvement was realized by incorporating wavelength linearization. Thus, for the low resolution first order flame




Figure 6 Plot of the absolute difference of the observed values and the fitted values versus wavelength for the conditions of Fig. 3.




Figure 7 Plot of the absolute difference of the observed and the fitted values versus wavelength for the conditions of Fig. 5.

data of Fig. 4, nonlinearity in the wavelength does not notably influence the quality of the fit.

Wavelength linearization has also been applied to the same kind of OH absorption data taken in propellant flames burning in 1.5 MPa nitrogen. This is a significantly more hostile environment than a premixed laminar flame, and there is a limited amount of time for accumulating the data. One of the best propellant spectra that we have been able to obtain is shown on Fig. 8. Although the second order propellant data is noisier than the flame data, the linearization is shown to still lower the uncertainty in the fitted values (see Table 2).




Figure 8 Rotationally resolved OH spectrum: 0.4 to 0.9 mm above the propellant surface, burning in 1.5 MPa nitrogen: data acquisition time, 1.6 sec.; monochromator, second order; data fitted using linearization.

Table 2. Conditions for which OH absorption spectra were least squares fitted to obtain temperatures and concentrations together with their standard deviation.

<u>Sample</u>	<u>Order</u>	<u>Wavelength</u>	<u>OH Concentration</u> (molecules/cc)	<u>Temperature</u> (K)	<u>Comments</u>
		<u>Linearization</u>			
Flame	2	no	$3.67 \pm 0.14 \times 10^{16}$	2214 ± 85	data and fit on Fig. 3
Flame	1	no	$3.05 \pm 0.08 \times 10^{16}$	2378 ± 64	data and fit on Fig. 4
Flame	2	yes	$3.66 \pm 0.07 \times 10^{16}$	2212 ± 37	shown on Fig. 5
Flame	1	yes	$3.07 \pm 0.07 \times 10^{16}$	2380 ± 64	same data as Fig. 3
					fit not shown
Propellant	2	no	$2.29 \pm 0.07 \times 10^{16}$	2346 ± 75	not shown
Propellant	2	yes	$2.33 \pm 0.06 \times 10^{16}$	2420 ± 54	data and fit on Fig. 8

V. CONCLUSION

The wavelength linearization process described in the text has significantly improved the fitting of data taken by a photodiode array. In the case of OH absorption spectra taken in a steady-state premixed flame, the standard deviation of the fit was lowered by more than a factor of two. As the spectral resolution decreases or other sources of noise become larger, the relative importance of the linearization effect decreases.

INTENTIONALLY LEFT BLANK.

REFERENCES

1. J.M. Katzenberger, Fran Adar, and J.M. Lerner, "An Improved Algorithm for Linearizing in Wavelength or Wavenumber Spectral Data Acquired with a Diode Array," *Proceedings of the Microbeam Analysis Conference*, p.165, San Francisco Press, 1987.
2. J. A. Vanderhoff and A. J. Kotlar, "Simultaneous Determination of Temperatures and OH Concentrations in a Solid Propellant Flame," *23rd Symposium (international) on Combustion*, to be published.
3. J. A. Vanderhoff and A. J. Kotlar, "Temperature and OH Concentrations in a Solid Propellant Flame Using Absorption Techniques," BRL Technical Report, BRL-TR-3098, Aberdeen Proving Ground, MD, April 1990
4. J. A. Vanderhoff, "Spectral Studies of Solid Propellant Combustion II. Emission and Absorption Results for M-30 and HMX1 Propellants", BRL Technical Report, BRL-TR-3055, Aberdeen Proving Ground, MD, December 1989.
5. J. A. Vanderhoff, "Species Profiles in Solid Propellant Flames Using Absorption and Emission Profiles," *Combust. Flame*, to be published.
6. J. M. Lerner and A. Thevenson, *The Optics of Spectroscopy: A Tutorial V2.0*, Instruments SA, Inc., 6 Olsen Avenue, Metuchen, NJ 08820, and Jobin-Yvon Division of Instruments SA, 16-18 rue du Canal, F-91160 Longjumeau, France.
7. R. P. Lucht, R. C. Peterson, and N. M. Laurendeau, "Fundamentals of Absorption Spectroscopy for Selected Diatomic Flame Radicals," PURDU-CL-78-06, 1978.
8. W. E. Wentworth, "Rigorous Least Squares Adjustment," *J. Chem. Ed.*, **42**, 96 (1965).

INTENTIONALLY LEFT BLANK.

<u>No. of</u> <u>Copies</u>	<u>Organization</u>	<u>No. of</u> <u>Copies</u>	<u>Organization</u>
1	Office of the Secretary of Defense OUSD(A) Director, Live Fire Testing ATTN: James F. O'Bryon Washington, DC 20301-3110	1	Director US Army Aviation Research and Technology Activity ATTN: SAVRT-R (Library) M/S 219-3 Ames Research Center Moffett Field, CA 94035-1000
2	Administrator Defense Technical Info Center ATTN: DTIC-DDA Cameron Station Alexandria, VA 22304-6145	1	Commander US Army Missile Command ATTN: AMSMI-RD-CS-R (DOC) Redstone Arsenal, AL 35898-5010
1	HQDA (SARD-TR) WASH DC 20310-0001	1	Commander US Army Tank-Automotive Command ATTN: AMSTA-TSL (Technical Library) Warren, MI 48397-5000
1	Commander US Army Materiel Command ATTN: AMCDRA-ST 5001 Eisenhower Avenue Alexandria, VA 22333-0001	1	Director US Army TRADOC Analysis Command ATTN: ATAA-SL White Sands Missile Range, NM 88002-5502
1	Commander US Army Laboratory Command ATTN: AMSLC-DL Adelphi, MD 20783-1145	(Class. only) 1	Commandant US Army Infantry School ATTN: ATSH-CD (Security Mgr.) Fort Benning, GA 31905-5660
2	Commander US Army, ARDEC ATTN: SMCAR-IMI-I Picatinny Arsenal, NJ 07806-5000	(Unclass. only) 1	Commandant US Army Infantry School ATTN: ATSH-CD-CSO-OR Fort Benning, GA 31905-5660
2	Commander US Army, ARDEC ATTN: SMCAR-TDC Picatinny Arsenal, NJ 07806-5000	1	Air Force Armament Laboratory ATTN: AFATL/DLODL Eglin AFB, FL 32542-5000
1	Director Benet Weapons Laboratory US Army, ARDEC ATTN: SMCAR-CCB-TL Watervliet, NY 12189-4050		<u>Aberdeen Proving Ground</u>
1	Commander US Army Armament, Munitions and Chemical Command ATTN: SMCAR-ESP-L Rock Island, IL 61299-5000	2	Dir, USAMSAA ATTN: AMXSY-D AMXSY-MP, H. Cohen
1	Commander US Army Aviation Systems Command ATTN: AMSAV-DACL 4300 Goodfellow Blvd. St. Louis, MO 63120-1798	1	Cdr, USATECOM ATTN: AMSTE-TO-F
		3	Cdr, CRDEC, AMCCOM ATTN: SMCCR-RSP-A SMCCR-MU SMCCR-MSI
		1	Dir, VIAMO ATTN: AMSLC-VL-D

<u>No. of</u> <u>Copies</u>	<u>Organization</u>
4	Commander US Army Research Office ATTN: R. Ghirardelli D. Mann R. Singleton R. Shaw P.O. Box 12211 Research Triangle Park, NC 27709-2211
2	Commander US Army, ARDEC ATTN: SMCAR-AEE-B, D.S. Downs SMCAR-AEE, J.A. Lannon Picatinny Arsenal, NJ 07806-5000
1	Commander US Army, ARDEC ATTN: SMCAR-AEE-BR, L. Harris Picatinny Arsenal, NJ 07806-5000
2	Commander US Army Missile Command ATTN: AMSMI-RK, D.J. Ifshin W. Wharton Redstone Arsenal, AL 35898
1	Commander US Army Missile Command ATTN: AMSMI-RKA, A.R. Maykut Redstone Arsenal, AL 35898-5249
1	Office of Naval Research Department of the Navy ATTN: R.S. Miller, Code 432 800 N. Quincy Street Arlington, VA 22217
1	Commander Naval Air Systems Command ATTN: J. Ramnarace, AIR-54111C Washington, DC 20360
1	Commander Naval Surface Warfare Center ATTN: J.L. East, Jr., G-23 Dahlgren, VA 22448-5000
2	Commander Naval Surface Warfare Center ATTN: R. Bernecker, R-13 G.B. Wilmot, R-16 Silver Spring, MD 20903-5000

<u>No. of</u> <u>Copies</u>	<u>Organization</u>
5	Commander Naval Research Laboratory ATTN: M.C. Lin J. McDonald E. Oran J. Shnur R.J. Doyle, Code 6110 Washington, DC 20375
1	Commanding Officer Naval Underwater Systems Center Weapons Dept. ATTN: R.S. Lazar/Code 36301 Newport, RI 02840
2	Commander Naval Weapons Center ATTN: T. Boggs, Code 388 T. Parr, Code 3895 China Lake, CA 93555-6001
1	Superintendent Naval Postgraduate School Dept. of Aeronautics ATTN: D.W. Netzer Monterey, CA 93940
3	AL/LSCF ATTN: R. Corley R. Geisler J. Levine Edwards AFB, CA 93523-5000
1	AI/MKPB ATTN: B. Goshgarian Edwards AFB, CA 93523-5000
1	AFOSR ATTN: J.M. Tishkoff Bolling Air Force Base Washington, DC 20332
1	OSD/SDIO/UST ATTN: L. Caveny Pentagon Washington, DC 20301-7100
1	Commandant USAFAS ATTN: ATSF-TSM-CN Fort Sill, OK 73503-5600
1	F.J. Seiler ATTN: S.A. Shackleford USAF Academy, CO 80840-6528

<u>No. of Copies</u>	<u>Organization</u>
1	University of Dayton Research Institute ATTN: D. Campbell AL/PAP Edwards AFB, CA 93523
1	NASA Langley Research Center Langley Station ATTN: G.B. Northam/MS 168 Hampton, VA 23365
4	National Bureau of Standards ATTN: J. Hastie M. Jacox T. Kashiwagi H. Semerjian US Department of Commerce Washington, DC 20234
1	Aerojet Solid Propulsion Co. ATTN: P. Micheli Sacramento, GA 95813
1	Applied Combustion Technology, Inc. ATTN: A.M. Varney P.O. Box 17885 Orlando, FL 32860
2	Applied Mechanics Reviews The American Society of Mechanical Engineers ATTN: R.E. White A.B. Wenzel 345 E. 47th Street New York, NY 10017
1	Atlantic Research Corp. ATTN: M.K. King 5390 Cherokee Avenue Alexandria, VA 22314
1	Atlantic Research Corp. ATTN: R.H.W. Waesche 7511 Wellington Road Gainesville, VA 22065
1	AVCO Everett Research Laboratory Division ATTN: D. Stickler 2385 Revere Beach Parkway Everett, MA 02149

<u>No. of Copies</u>	<u>Organization</u>
1	Battelle Memorial Institute Tactical Technology Center ATTN: J. Huggins 505 King Avenue Columbus, OH 43201
1	Cohen Professional Services ATTN: N.S. Cohen 141 Channing Street Redlands, CA 92373
1	Exxon Research & Eng. Co. ATTN: A. Dean Route 22E Annandale, NJ 08801
1	Ford Aerospace and Communications Corp. DIVAD Division Div. Hq., Irvine ATTN: D. Williams Main Street & Ford Road Newport Beach, CA 92663
1	General Applied Science Laboratories, Inc. 77 Raynor Avenue Ronkonkoma, NY 11779-6649
1	General Electric Armament & Electrical Systems ATTN: M.J. Bulman Lakeside Avenue Burlington, VT 05401
1	General Electric Ordnance Systems ATTN: J. Mandzy 100 Plastics Avenue Pittsfield, MA 01203
2	General Motors Rsch Labs Physics Department ATTN: T. Sloan R. Teets Warren, MI 48090
2	Hercules, Inc. Allegheny Ballistics Lab. ATTN: W.B. Walkup E.A. Yount P.O. Box 210 Rocket Center, WV 26726

<u>No. of</u> <u>Copies</u>	<u>Organization</u>
1	Honeywell, Inc. Government and Aerospace Products ATTN: D.E. Broden/ MS MN50-2000 600 2nd Street NE Hopkins, MN 55343
1	Honeywell, Inc. ATTN: R.E. Tompkins MN38-3300 10400 Yellow Circle Drive Minnetonka, MN 55343
1	IBM Corporation ATTN: A.C. Tam Research Division 5600 Cottle Road San Jose, CA 95193
1	IIT Research Institute ATTN: R.F. Remaly 10 West 35th Street Chicago, IL 60616
2	Director Lawrence Livermore National Laboratory ATTN: C. Westbrook M. Costantino P.O. Box 808 Livermore, CA 94550
1	Lockheed Missiles & Space Co. ATTN: George Lo 3251 Hanover Street Dept. 52-35/B204/2 Palo Alto, CA 94304
1	Los Alamos National Lab ATTN: B. Nichols T7, MS-B284 P.O. Box 1663 Los Alamos, NM 87545
1	National Science Foundation ATTN: A.B. Harvey Washington, DC 20550
1	Olin Corporation Smokeless Powder Operations ATTN: V. McDonald P.O. Box 222 St. Marks, FL 32355

<u>No. of</u> <u>Copies</u>	<u>Organization</u>
1	Paul Gough Associates, Inc. ATTN: P.S. Gough 1048 South Street Portsmouth, NH 03801-5423
2	Princeton Combustion Research Laboratories, Inc. ATTN: M. Summerfield N.A. Messina 475 US Highway One Monmouth Junction, NJ 08852
1	Hughes Aircraft Company ATTN: T.E. Ward 8433 Fallbrook Avenue Canoga Park, CA 91303
1	Rockwell International Corp. Rocketdyne Division ATTN: J.E. Flanagan/HB02 6633 Canoga Avenue Canoga Park, CA 91304
4	Sandia National Laboratories Division 8354 ATTN: R. Cattolica S. Johnston P. Mattern D. Stephenson Livermore, CA 94550
1	Science Applications, Inc. ATTN: R.B. Edelman 23146 Cumorah Crest Woodland Hills, CA 91364
3	SRI International ATTN: G. Smith D. Crosley D. Golden 333 Ravenswood Avenue Menlo Park, CA 94025
1	Stevens Institute of Tech. Davidson Laboratory ATTN: R. McAlevy, III Hoboken, NJ 07030
1	Thiokol Corporation Elkton Division ATTN: S.F. Palopoli P.O. Box 241 Elkton, MD 21921

No. of
Copies Organization

- 1 Morton Thiokol, Inc.
Huntsville Division
ATTN: J. Deur
Huntsville, AL 35807 7501
- 3 Thiokol Corporation
Wasatch Division
ATTN: S.J. Bennett
P.O. Box 524
Brigham City, UT 84302
- 1 United Technologies Research Center
ATTN: A.C. Eckbreth
East Hartford, CT 06108
- 3 United Technologies Corp.
Chemical Systems Division
ATTN: R.S. Brown
T.D. Myers (2 copies)
P.O. Box 49028
San Jose, CA 95161-9028
- 1 Universal Propulsion Company
ATTN: H.J. McSpadden
Black Canyon Stage 1
Box 1140
Phoenix, AZ 85029
- 1 Veritay Technology, Inc.
ATTN: E.B. Fisher
4845 Millersport Highway
P.O. Box 305
East Amherst, NY 14051-0305
- 1 Brigham Young University
Dept. of Chemical Engineering
ATTN: M.W. Beckstead
Provo, UT 84058
- 1 California Institute of Tech.
Jet Propulsion Laboratory
ATTN: L. Strand/MS 512/102
4800 Oak Grove Drive
Pasadena, CA 91009
- 1 California Institute of
Technology
ATTN: F.E.C. Culick/
MC 301-46
204 Karman Lab.
Pasadena, CA 91125
- 1 University of California
Los Alamos Scientific Lab.
P.O. Box 1663, Mail Stop B216
Los Alamos, NM 87545

No. of
Copies Organization

- 1 University of California,
San Diego
ATTN: F.A. Williams
AMES, B010
La Jolla, CA 92093
- 2 University of California,
Santa Barbara
Quantum Institute
ATTN: K. Schofield
M. Steinberg
Santa Barbara, CA 93106
- 1 University of Colorado at
Boulder
Engineering Center
ATTN: J. Daily
Campus Box 427
Boulder, CO 80309-0427
- 2 University of Southern
California
Dept. of Chemistry
ATTN: S. Benson
C. Wittig
Los Angeles, CA 90007
- 1 Case Western Reserve Univ.
Div. of Aerospace Sciences
ATTN: J. Tien
Cleveland, OH 44135
- 1 Cornell University
Department of Chemistry
ATTN: T.A. Cool
Baker Laboratory
Ithaca, NY 14853
- 1 University of Delaware
ATTN: T. Brill
Chemistry Department
Newark, DE 19711
- 1 University of Florida
Dept. of Chemistry
ATTN: J. Winetfordner
Gainesville, FL 32611
- 3 Georgia Institute of
Technology
School of Aerospace
Engineering
ATTN: E. Price
W.C. Strahle
B.T. Zinn
Atlanta, GA 30332

<u>No. of</u> <u>Copies</u>	<u>Organization</u>
1	University of Illinois Dept. of Mech. Eng. ATTN: H. Krier 144MEB, 1206 W. Green St. Urbana, IL 61801
1	Johns Hopkins University/API Chemical Propulsion Information Agency ATTN: T.W. Christian Johns Hopkins Road Laurel, MD 20707
1	University of Michigan Gas Dynamics Lab Aerospace Engineering Bldg. ATTN: G.M. Faeth Ann Arbor, MI 48109-2140
1	University of Minnesota Dept. of Mechanical Engineering ATTN: E. Fletcher Minneapolis, MN 55455
3	Pennsylvania State University Applied Research Laboratory ATTN: K.K. Kuo H. Palmer M. Micci University Park, PA 16802
1	Pennsylvania State University Dept. of Mechanical Engineering ATTN: V. Yang University Park, PA 16802
1	Polytechnic Institute of NY Graduate Center ATTN: S. Lederman Route 110 Farmingdale, NY 11735
2	Princeton University Forrestal Campus Library ATTN: K. Brezinsky I. Glassman P.O. Box 710 Princeton, NJ 08540
1	Purdue University School of Aeronautics and Astronautics ATTN: J.R. Osborn Grissom Hall West Lafayette, IN 47906

<u>No. of</u> <u>Copies</u>	<u>Organization</u>
1	Purdue University Department of Chemistry ATTN: E. Grant West Lafayette, IN 47906
2	Purdue University School of Mechanical Engineering ATTN: N.M. Laurendeau S.N.B. Murthy TSPC Chaffee Hall West Lafayette, IN 47906
1	Rensselaer Polytechnic Inst. Dept. of Chemical Engineering ATTN: A. Fontijn Troy, NY 12181
1	Stanford University Dept. of Mechanical Engineering ATTN: R. Hanson Stanford, CA 94305
1	University of Texas Dept. of Chemistry ATTN: W. Gardiner Austin, TX 78712
1	University of Utah Dept. of Chemical Engineering ATTN: G. Flandro Salt Lake City, UT 84112
1	Virginia Polytechnic Institute and State University ATTN: J.A. Schetz Blacksburg, VA 24061
1	Freedman Associates ATTN: E. Freedman 2411 Diana Road Baltimore, MD 21209-1525

USER EVALUATION SHEET/CHANGE OF ADDRESS

This Laboratory undertakes a continuing effort to improve the quality of the reports it publishes. Your comments/answers to the items/questions below will aid us in our efforts.

1. BRL Report Number BRL-MR-3866 Date of Report SEPTEMBER 1990
2. Date Report Received _____
3. Does this report satisfy a need? (Comment on purpose, related project, or other area of interest for which the report will be used.) _____

4. Specifically, how is the report being used? (Information source, design data, procedure, source of ideas, etc.) _____

5. Has the information in this report led to any quantitative savings as far as man-hours or dollars saved, operating costs avoided, or efficiencies achieved, etc? If so, please elaborate. _____

6. General Comments. What do you think should be changed to improve future reports? (Indicate changes to organization, technical content, format, etc.) _____

CURRENT
ADDRESS

Name

Organization

Address

City, State, Zip Code

7. If indicating a Change of Address or Address Correction, please provide the New or Correct Address in Block 6 above and the Old or Incorrect address below.

OLD
ADDRESS

Name

Organization

Address

City, State, Zip Code

(Remove this sheet, fold as indicated, staple or tape closed, and mail.)

-----FOLD HERE-----

DEPARTMENT OF THE ARMY

Director
U.S. Army Ballistic Research Laboratory
ATTN: SLCBR-DD-T
Aberdeen Proving Ground, MD 21005-9989
OFFICIAL BUSINESS



NO POSTAGE
NECESSARY
IF MAILED
IN THE
UNITED STATES

BUSINESS REPLY MAIL
FIRST CLASS PERMIT No 0001, APG, MD

POSTAGE WILL BE PAID BY ADDRESSEE

Director
U.S. Army Ballistic Research Laboratory
ATTN: SLCBR-DD-T
Aberdeen Proving Ground, MD 21005-9989



-----FOLD HERE-----